

DLD 8080



Custom sized Delayline Detector with enlarged active area

Features



- CF150 mounting flange
 - 80 x 80 mm² active area of DLD body and \varnothing 82 mm active MCP area
 - Linear response due to single event counting
 - Extremely low dark count rate: < 20 cps
 - Up to 3 MHz count rate in 2D mode
 - Up to 41 μ m of pixel size
 - Up to 81 ps digital time bin resolution
 - Time histogramming for time of flight (tof) measurements
 - Also available as hybrid version with an additional integrated 9 channel detector
-
- For further information contact Surface Concept GmbH